

SERASYNTH

Polydioxanone monofilament, violet and undyed.

The material is marked not only by its outstanding tissue glide, but also its linear and knot tensile strength. Serasynth possesses prolonged absorption behaviour.

- For indications requiring a long-term wound healing phase

Approximate Absorption times

- Half-value 28-42 d, dissolution time 180-210 d

Benefits

- Excellent sliding capabilities
- No capillarity
- Reducing the risk of infection
- Optimal handling due to approved Alu-peel packaging

Indications

- Fascia closure
- Abdominal wall closure
- Stomach and small intestine
- Large intestine
- Kidneys
- Urinary tract
- Thorax closure
- Tendons
- Ligaments
- Semilunar cartilage
- Bone



*Ethical Agents Veterinary Marketing Ltd, 54 Hobill Ave, Manukau, 2104, New Zealand.
P O Box 97110, Manukau City, Manukau, 2241, New Zealand.
Telephone 64-9-262 1388, Fax 64-9-262 1411, Freephone 0800 800 624.
E-mail eage@xtra.co.nz, Website www.ethicalagents.co.nz.*